

**SPECIFICATIONS**

| Part No. | Alloy | Bolt Type | Product Weight | Max Patient Weight | Screw Torque (M8x12mm) |
|----------|----------|-----------|----------------|--------------------|------------------------|
| 2T09 | Titanium | Steel M8 | 80 gr. | 100 kg. / 220 lbs. | 15 Nm. |
| 2T10 | Titanium | Steel M10 | 93 gr. | 136 kg. / 300 lbs. | 15 Nm. |
| 2T11 | Titanium | Steel M10 | 80 gr. | 136 kg. / 300 lbs. | 15 Nm. |
| 2S10 | S. Steel | Steel M10 | 130 gr. | 136 kg. / 300 lbs. | 15 Nm. |
| 2A10 | Aluminum | Steel M10 | 106 gr. | 100 kg. / 220 lbs. | 15 Nm. |

USAGE

This product is designed to be used for lower limb external prosthetic applications. Before fitting this product onto the patient, all parts must be joined properly and all set screws must be secured enough for the patient to walk. In order to make final adjustments, there should be no loose screws.

ATTENTION

After the final adjustments, apply single drop of removable Loctite® 242 Blue® to each of the set/cap screws and then torque the product. Improper installation or uneven, insufficient or excessive torque may cause bolt or adapter failure. Always use quality metric tools for installation.

WARRANTY

All MEDEX Modular Adapters are unconditionally guaranteed for 2 years. Cutting, drilling, modifying or using the component other than as described in this manual will void the warranty.