

**SPECIFICATIONS**

Titanium, stainless steel and aluminum tube clamp adapters for 30 mm pylon tubes.

Part No.	Alloy	Product Weight	Max Patient Weight	Set Screw Torque	Clamp Screw Torque
2T20	Titanium	70 gr.	136 Kg. / 300 lbs.	15 Nm.	9 Nm.
2S20	S. Steel	113 gr.	136 Kg. / 300 lbs.	15 Nm.	9 Nm.
2A20	Aluminum	68 gr.	100 Kg. / 220 lbs.	15 Nm.	9 Nm.

USAGE

This product is designed to be used for lower limb external prosthetic toe direction and pylon height adjustment applications. Before fitting this product onto the patient, all parts must be joined properly and all set screws must be secured enough for the patient to walk. The pylon tube has to be pushed into this product until it sits on the adapter rim. In order to make final adjustments, there should be no loose screws.

ATTENTION

After the final adjustments, apply single drops of removable Loctite® [222R] to each of the set/cap screws and then torque the product. Improper installation or uneven, insufficient or excessive torque may cause bolt or adapter failure. Always use quality metric tools for installation.

WARRANTY

All MEDEX Modular Adapters are unconditionally guaranteed for 2 years. Cutting, drilling, modifying or using the component other than as described in this manual will void the warranty.