**Technical Specification**

Part No.	Alloy	Product Weight	Max Patient Weight
5T40	Titanium	349 gr.	100 Kg. / 220 lbs.
5S40	Aluminum	525 gr.	100 Kg. / 220 lbs.

**Description and Purpose**

These instructions are for use by the practitioner. This product is designed to be used for lower limb external prosthetic length, toe direction and alignment corrections. Intended for single user only. This product can be used for mobility classes 1-4. There are no known contra-indications if used in accordance with these instructions. Ensure that the user has understood any instructions for use, drawing particular attention to the maintenance information.

**Usage**

The large screw on the swing block is structured to work in conjunction with pressure plate's adjustment screw which is mounted on plastic, under the proximal pyramid. It is not advisable to tighten the swing block screw while the pressure plate screw is loose and vice versa.

In order to make adjustments to the knee joint:

1. If brake is too soft and cannot achieve proper movement on the knee joint, start by turning the swing block screw counter-clockwise until it becomes non-effective.
2. Begin rotating the pressure plate screw clockwise until the knee joint locks in. At this moment the knee joint will not move freely.
3. Turn pressure plate screw counter clockwise until the knee joints begins moving with resistance. Once this stiff movement is achieved, stop rotating the pressure plate screw.
4. Rotate the swing block screw clockwise until there is slightly less resistance in the movement of the knee joint.

There will be noticeably less friction at this point.

5. Rotate the pressure plate screw counter clockwise until the desired brake safety and friction is achieved.

For very active patients or moderate activity patients, the brake catch should not be very sensitive. Turning the swing block screw clockwise ¼ or ½ turn should be sufficient for this activity level. Always verify that the brake is in working condition and will engage. For patients with low activity level, this final adjustment is not necessary.

**Attention**

Final adjustments should be made within 3-5 days.

**Limitations on Use**

Intended life:

A local risk assessment should be carried out based upon activity and usage.

Lifting loads:

User weight and activity is governed by the stated limits. Load carrying by the user should be based on a local risk assessment.

Environment:

Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts. Thoroughly rinse with fresh water after use in salt or chlorinated water. Exclusively for use between -15°C and 50°C (5°F to 122°F).

**Maintenance**

Maintenance is to be carried out by competent personnel only. An annual visual inspection is recommended. During maintenance, it is appropriate to consider any weight or activity increase by the user. If the prosthesis is exposed to circumstances or impact for which it was not designed, inspection by a qualified practitioner is recommended.

- Ensure all screws are secure. If not, remove and clean screws, reapply Loctite to the pyramid screws only, tighten all screws to the correct torque settings.
- Check for visual defects that may affect proper function.

The user should be advised to contact their practitioner if:

- Their condition changes.
- Performance of the device changes.
- Changes in performance may include:
  - Any unusual noises
  - Excessive play or loss of alignment

**Warranty**

All MEDEX Modular Adapters are unconditionally guaranteed for 2 years. Cutting, drilling, modifying or using the component other than as described in this manual will void the warranty.

**Liability**

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device.

**CE Conformity**

This product meets the requirements of MDR 2017/745 guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines.

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Teile-Nr.	Legierung	Produktgewicht	Maximales Patientengewicht
5T40	Titan	349 gr.	100 Kg. / 220 lbs.
5S40	Aluminium	525 gr.	100 Kg. / 220 lbs.

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