Technical Specification

Part No.	Alloy	Product Weight	Build Height	Max Patient Weight	Torque
4A20	Titanium / Aluminum / S. Steel	130 g	47 mm.	136 Kg. / 300 lbs.	9 Nm.
4T20	Titanium / Aluminum / S. Steel	139 g	49 mm.	136 Kg. / 300 lbs.	9 Nm.
4T21	Titanium / Aluminum / S. Steel	146 g	49 mm.	136 Kg. / 300 lbs	9 Nm.

USAGE

Cast Preparation and RIM Installation:

1. Use RIM as the lamination dummy. Remove lock from the RIM. Tape both sides of the RIM to block resin flow inside. Install dummy screws (3) to block resin. If the dummy screws are not available, use silicon putty to protect the threaded holes from resin. After the final cast alteration, place the RIM on the cast and check for close and proper fit with RIM at the distal end of the socket. If RIM does not fit to cast properly, make alterations on the cast to achieve a proper form. Do not continue to next section until best possible fit is achieved.
2. After determining RIM alignment based on patient characteristics, roughen the cast at the distal end where RIM will be applied. Install RIM using plaster of Paris after desired alignment. Scrape away any excess plaster. Fill up joint area to assure smooth interface between cast and RIM. Once the plaster has dried, sand the cast around the RIM. Then twist off and remove the RIM from the cast. Slice off any excess plaster to leave a smooth rounded surface. With the RIM removed, pull the PVA bag tightly over the cast. Use a heat gun to shrink the bag around the cast ensuring that no bumps or air bubbles remain. After attaching the proximal end of the bag to vacuum suction, twist and clamp bag at the distal end of the cast. Cut off any excess above the clamp. Tape the proximal end of the bag tightly below the cast to ensure the bag is fully sealed. Remove the clamp from distal end and apply tape. This will help protect the inner bag from punctures during the lamination process. Place the RIM back on the cast. Use rolled silicon to seal the edges of the RIM to ensure no lamination fluid can flow between the RIM and the tape. Pull nylon stockinet over the RIM and the cast. Tie a string around the RIM groove, and then pull the stockinet down and over the cast.
3. Next, pull a carbon fiber braid over the cast, tie a string on the RIM groove, and pull the braid down and over the cast. Build enforcement layers as necessary. Add another nylon stockinet over the carbon fiber braid, tie a string around the RIM groove and pull back down over the cast. Now add a PVA bag attaching the lower end of the bag to suction. Use a heat gun to shrink the bag around the cast ensuring that no bumps or air bubbles remain. Mix the resin and pour it inside the PVA from the top and squeeze the bag to push the resin on to the cast. Smooth the cast using a string to ensure uniform layer of laminate across the surface.

Lock Assembly:

1. Remove socket from cast, clean up edges, unscrew lamination bolts and remove tape and silicon from threaded holes of the RIM and install (3) M8 socket set screws. Insert shuttle lock through the hole and secure rotation with the M8 screws. Do not torque. Please make sure that the entire assembly area is free of dirt or debris that may stop the lock from rotating.
2. After installing shuttle lock unit into its location, check for proper fit. Two pieces should now be in full contact with no space between. Check from inside the socket as well to ensure smooth engagement between the RIM and lock.
3. Screw on rotation control screws until they engage with the inner core. Do not torque. Turn them counter clockwise one quarter turn to disengage with inner core. The lock is now ready for rotation adjustment.
4. Rotation stop with one screw during fitting is optional, however you must torque all rotation set screws to 9 Nm. after the final adjustment and before use by patient.

Maintenance:

Maintenance is to be carried out by competent personnel only. An annual visual inspection is recommended. During maintenance, it is appropriate to consider any weight or activity increase by the user. If the prosthesis is exposed to circumstances or impact for which it was not designed, inspection by a qualified practitioner is recommended.

- Ensure all screws are secure. If not, remove and clean screws, reapply Loctite to the pyramid screws only, tighten all screws to the correct torque settings.
- Check for visual defects that may affect proper function. The user should be advised to contact their practitioner if:
 - Their condition changes.
 - Performance of the device changes.
 - Changes in performance may include:
 - o Any unusual noises
 - o Excessive play or loss of alignment

Warranty:

All MEDEX Modular Adapters are unconditionally guaranteed for 2 years. Cutting, drilling, modifying or using the component other than as described in this manual will void the warranty.

Liability:

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device.

CE Conformity:

This product meets the requirements of MDR 2017/745 guidelines for medical products. This product has been classified as a Class 1 Product according to the classification criteria outlined in Appendix IX of the guidelines.

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